

Single-Use Bayonet Style Doppler Probes

Indicated for use to obtain blood flow sounds from intraoperatively exposed vessels of interest, thereby verifying patency or revealing absent or significantly diminished flow sounds.

- Pre-sterile and packaged for single use
- No sterilization processing time; ready to use
- Eliminates the risk of cross contamination
- Easy and convenient for busy hospitals

MODELS AVAILABLE

NDP-10 10 MHz Single-Use Bayonet Style Doppler Probe (2.5mm diameter)
NDP-20 20 MHz Single-Use Bayonet Style Doppler Probe (2.5mm diameter)

USE WITH HADECO DOPPLERS



ES-100VX MINIDOPPLER®



BIDOP® 3



SMARTDOP® 45

Single-use probes are sterile and nonpyrogenic (unless package is opened or damaged). Not for Implantation. See the operation manual provided with your Doppler for any Doppler cautions & warnings.

MiniDoppler, Bidop, and Smartdop are registered trademarks of Koven Technology, Inc., St. Louis, MO USA

VASCULAR INSTRUMENTS
KOVEN
TECHNOLOGY INCORPORATED
Your Vascular Healthcare Partner

KOVEN TECHNOLOGY, INC.
477 N. LINDBERGH BLVD., SUITE 220
ST. LOUIS, MO 63141
T: +1.314.542.2101
WWW.KOVEN.COM

UK
CA 0086

CE 2797

DS-1660 rev. 1

Made in Japan/Packaged, Sterilized in USA