

Single-Use Vascular Doppler Probes

For use with Hadeco Dopplers

Indicated for use to obtain blood flow sounds from intraoperatively exposed vessels of interest, thereby verifying patency or revealing absent or significantly diminished flow sounds.

- Pre-sterile and packaged for single use
- No sterilization processing time; ready to use
- Eliminates the risk of cross contamination

Bayonet Style Single-Use Probes

Pencil Style Single-Use Probes

MODELS AVAILABLE

VDP-8 8 MHz Single-Use Pencil Style Doppler Probe (4.8mm diameter)

NDP-10 10 MHz Single-Use Bayonet Style Doppler Probe (2.5mm diameter)

NDP-20 20 MHz Single-Use Bayonet Style Doppler Probe (2.5mm diameter)

*Single-use probes are sterile and nonpyrogenic (unless package is opened or damaged).
Not for Implantation.

See the operation manual provided with your Doppler for any Doppler cautions & warnings.