

# VDP-8

## Single Use Vascular Doppler Probes

Indicated for use to obtain blood flow sounds from intraoperatively exposed vessels of interest, thereby verifying patency or revealing absent or significantly diminished flow sounds.

- Pre-sterile and packaged for single use
- No sterilization processing time; ready to use
- Eliminates the risk of cross contamination
- Easy and convenient for busy hospitals



USE WITH  
HADECO VASCULAR DOPPLERS 



ES-100VX MINIDOPPLER



BIDOP 3

See the operation manual provided with your Doppler for any Doppler cautions & warnings.

VDP-8 8 MHz probe diameter 4.8mm

Single-use probes are sterile and nonpyrogenic (unless package is opened or damaged). Not for Implantation.

VASCULAR INSTRUMENTS  
**KOVEN**  
TECHNOLOGY INCORPORATED  
Your Vascular Healthcare Partner

**KOVEN TECHNOLOGY, INC.**  
477 N. LINDBERGH BLVD., SUITE 220  
ST. LOUIS, MO 63141  
T: +1.314.542.2101  
WWW.KOVEN.COM

**UK**  
**CA** 0086  2797

DS-1661 rev. 1

Made in Japan/Packaged, Sterilized in USA