


SAUGERMED


Surgical Brain Strips

CE₁₇₈₃

7 different strips for different purposes.



- Produced by nonwoven form with materials of hydrofiber speciality.
- Atraumatic and soft.
- Unfeathered, unfuzzy, does not leave loofah.
- They can haemostatize smallest compressions.
- All strips contains x-ray detectable Radio-pack speciality blue strap
- Strips are attached to the counter card, which consists a protective. band made of carton to prevent from damage.



1:1 scale product table (mm.)

*Please contact with our sales department for order informations of Saugermед.

Surgical Brain Strips



SAUGERMED Thinpad

- Absorbation capacity 8 times of its weight.
- Very thin, smooth and silky structured.
- Made of viscose and synthetic additive material.



SAUGERMED Plus

- Absorbation capacity 11 times of its weight .
- Thin, closely structured.
- Made of viscose and synthetic additive material.



SAUGERMED Thincot

- Absorbation capacity 9 times of its weight.
- Thin, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.



SAUGERMED Laminat

- Absorbation capacity 10 times of its weight.
- Standard thickness, smooth structured, laminated surface.
- Made of viscose and synthetic additive material.




SAUGERMED Soft

- Absorbation capacity 10 times of its weight
- Standard thickness, closely structured, laminated surface.
- Made of viscose and synthetic additive material.



SAUGERMED Ultra Soft

- Absorbation capacity 15 times of its weight.
- Smooth , delicate structured.
- Made of viscose and synthetic additive material.



SAUGERMED Cotton Soft

- Absorbation capacity 16 times of its weight.
- Very smooth, delicate structured.
- Made of viscose and synthetic additive material.

